

Plant protoplasts

PROTOPLASTS are cells (plant, fungal or bacterial) that have had their cell walls removed. This can be done mechanically, or by enzymic digestion. The 'naked' cells are surrounded only by a cell membrane and can be used in a variety of ways. For example, two or more protoplasts can be fused with the help of a detergent, polyethylene glycol, to produce hybrid cells with characteristics from each 'parent'. Infection of protoplasts with genetically-modified *Agrobacterium tumefaciens* is one way of introducing new genes into plant cells. Whole plants can be regenerated from protoplasts grown on solid or liquid media.

A mixture of carbohydrases can be used to degrade plant tissue, producing a protoplast suspension which can easily be seen under a microscope. For school work, this procedure does not require an aseptic environment if the protoplasts are not to be cultured.

Materials

Round green lettuce, 1 leaf
 21% sucrose solution, 0.1 cm³
 13% sorbitol solution, 20 cm³
 (sorbitol is available from high street chemists as an alternative to sugar for diabetics e.g. *Sionon*)
 Carbohydrase enzyme mixture, Novo *Viscozyme*[®], 0.5 cm³ (available from the NCBE)
 Paper tissue
 10 cm³ plastic syringe } for dispensing liquids
 1 cm³ plastic syringes, 2 }
 Test tube
 Centrifuge tube
 Fine glass stirring rod
 Filter funnel
 Strip of nylon gauze, to plug filter funnel
 Microscope slide
 Cover slip
 Microscope with x40 objective
 Water bath set at 37°C
 Bench centrifuge

Practical details

Lettuce preparation

1. Cut the lettuce leaf into pieces roughly 5 mm x 5 mm.
2. Add 15–20 lettuce pieces to 9.5 cm³ 13% sorbitol solution in a test tube.
3. Incubate the tube in a water bath at 37°C for 5 minutes. *During this time the lettuce tissue will equilibrate to the osmotic potential of the sorbitol solution.*

Enzymic digestion of cell walls

1. Gently stir 0.5 cm³ of the carbohydrase enzyme mixture into the sorbitol and lettuce preparation.
2. Return the tube to the water bath for another 20 minutes. Gently agitate the tube from time to time while it is incubating. *The carbohydrase mixture contains pectinases, cellulases and other cell wall-degrading enzymes.*

Recovery of protoplasts

1. Tightly pack the spout of the filter funnel with the nylon gauze. *This will produce a fine filter of roughly the correct pore size (able to trap cell debris greater than 50µm in diameter).*
2. Pour the digested lettuce suspension into the filter funnel.
3. Wash any trapped protoplasts through the filter using 10 cm³ of 13% sorbitol solution. Collect all the filtrates in a centrifuge tube. *The protoplasts will pass through the filter while the majority of the undigested cell debris will remain on top of the gauze.*
4. Ensure that the centrifuge is correctly balanced then spin the filtrate for roughly 5 minutes at 2,000 rpm (you may need to adjust this according to the type of centrifuge available).
5. Carefully discard the supernatant, leaving the pellet of protoplasts at the bottom of the tube. Use a tissue to wipe away excess liquid from the inside wall of the tube.
6. Resuspend the pellet in approximately 0.1 cm³ of 21% sucrose solution. *This is an osmotically-balanced (isotonic) medium in which the delicate protoplasts will not burst.*

Examination of protoplasts

The protoplasts will remain stable in the sucrose solution for several hours and can easily be seen without staining using a microscope with a x40 objective.

Safety

Care should be exercised when loading and operating the centrifuge. Students using the centrifuge must be closely supervised. Spills of enzyme should be wiped up promptly using plenty of water.


Further activities

Plant protoplasts are a safe alternative to blood cells for observing crenation and osmotically-induced lysis. Students can investigate the effects of a range of different concentrations of sucrose solutions on the appearance of the protoplasts.



ADDITIONAL INFORMATION

Plant protoplasts by Teresa Bengochea and John Dodds (1986) Chapman and Hall.
 ISBN: 0 412 26640 7.



Plant protoplasts

